

Mikrobac[®] Tissues

Alcohol-free disinfection wipes in the convenient flow pack. For the cleaning disinfection of alcohol-sensitive surfaces and medical devices.



Mikrobac[®] Tissues

Advantages at a glance

- ready-to-use disinfection wipes
- easy and safe dispensing of individual wipes out of resealable flow pack
- alcohol-, aldehyde, colourant and fragrance-free
- patented system solution: optimal release of active ingredients during disinfection
- also suitable for alcohol-sensitive materials
- can be used alternating with products based on other active substances (e.g. aldehydes or amines) without any problems

Application

Mikrobac Tissues are suitable for the convenient cleaning disinfection of sensitive surfaces and medical devices.

Areas of application

Mikrobac Tissues are suitable for the uncomplicated cleaning disinfection of:

Medical equipment as defined by the Medical Device Directive (MDD), such as:

- alcohol-sensitive application parts that cannot be submerged in liquids, such as ultrasonic heads of probes for abdominal examinations.

All kinds of alcohol sensitive medical inventory and surfaces, such as:

- surfaces in hospitals, medical and dental practices, ambulances
- surfaces in medical laboratories, residential and nursing homes
- small objects and surfaces such as toilet seats, door handles, bed frames, and office desks

Directions for use

Carefully wipe surfaces with Mikrobac Tissues. Ensure complete coverage so that optimum disinfection is guaranteed. If necessary, wipe off surface with a single use cloth after expiry of exposure time (e.g. if people are likely to come into direct contact with the surface). Dispose of tissue appropriately after use.

We recommend wearing suitable gloves while using this product. Not to be used for skin cleansing.

To remove disinfectant residue from sensitive plastic surfaces of medical devices, wipe surfaces with a cloth soaked in water (at least drinking water quality) after the exposure time.

When switching from one product to another, an intermediate cleaning has to be carried out.

Please follow the manufacturer's instructions. Do not allow disinfection solution to get inside of electrical devices. Not to be used for the final disinfection of semi-critical medical devices (e.g. ultrasonic heads of probes for vaginal examinations).

In order to prevent the wipes from drying out close lid when not in use.

Use disinfectants safely. Always read label and product information before use.



Mikrobac® Tissues

Proven efficacy

	Condition	Exposure time
Certified application recommendation (based on suspension and practical tests of VAH - German Association for Applied Hygiene)	Bactericidal / Yeasticidal dirty	5 min
Efficacy against enveloped viruses (DVV - German Society for the Control of Viral Diseases)	Virucidal against enveloped viruses (incl. HBV, HIV, HCV) (DVV / RKI)	30 sec

Bacteria and Fungi

Efficacy according to EN Phase 2 / Step 1 (suspension test)	Bactericidal (EN 13727)	dirty	30 sec
	Yeasticidal (EN 13624)	dirty	30 sec

Viruses

Efficacy against non-enveloped viruses	Norovirus, murine (EN 14476)	clean + dirty	4h
Efficacy against non-enveloped viruses (DVV - German Society for the Control of Viral Diseases)	Rotavirus (DVV / RKI)		30 sec
	Polyomavirus (DVV / RKI)		1 min

Food

Efficacy according to EN Phase 2 / Step 2 and Phase 2 / Step 1	Bactericidal / Yeasticidal (EN 13697 + EN 1276 + EN 1650)	low (20°C)	1 min
		low (10°C)	1 min
		low (4°C)	1 min

Microbiology

- bactericidal
- yeasticidal
- virucidal against enveloped viruses (incl. HBV, HIV, HCV)
- MNV
- rota- and polyomavirus

Composition

Active ingredients of the solution:
Benzyl-C12-18-alkyldimethyl-ammonium chloride
4 mg/g; didecyldimethylammonium chloride
4 mg/g

Chemical-physical data

All data referring to the wipes:

Dimensions:

Flow pack with 80 tissues: 180 x 200mm

Flow pack with 40 XXL tissues: 250 x 380mm

Reach:

Tissue 180 x 200 mm: approx. 0,85 m²

Tissue 250 x 380 mm: approx. 1,64 m²

The reach of Mikrobac Tissues not only depends on their size, but also on further factors, for example, the ambient temperature and structure of the surface to be disinfected. Always make sure to completely cover the surface area.

Related Products

Mikrobac forte+BODE X-Wipes

- **Mikrobac forte:** Aldehyde-free surface disinfectant cleaner with material-friendly protection factor.
- **BODE X-Wipes:** All-purpose fleece wipe dispenser system for all liquid HARTMANN surface disinfectants. Patented System solution: optimal active substance release during disinfection.
- **Bacillol 30 Tissues:** Material-friendly rapid disinfection wipes for sensitive surfaces in a convenient flow pack.



Product Presentation

Product	Content	Item no.
Mikrobac® Tissues		
	Flowpack with 80 tissues	on request
	Flowpack with 40 XXL tissues	on request

Please note: that the availability of products in the Mikrobac range may vary in different countries and regions. Contact your local distribution partner for more information.

The recommendations regarding our preparations are based on scientific tests and are given in good faith. More detailed recommendations, e.g. regarding material compatibility, are possible only in separate, individual cases. Our recommendations are not binding and do not constitute a guarantee. They do not preclude a company's own testing for the intended purpose and process. In this respect we cannot accept any liability. This is in accordance with our general conditions of sale and supply.

Manufacturer

BODE Chemie GmbH

A company of the HARTMANN GROUP

Melanchthonstraße 27

22525 Hamburg

Germany