

BOSCH

Bosch Healthcare Solutions GmbH
Stuttgarter Straße 130
71332 Waiblingen
Germany



CE
0123

IVD

F09G100186_09
Latest change: 04.05.2020



Bosch Healthcare Solutions GmbH
Stuttgarter Straße 130
71332 Waiblingen, Germany



BOSCH

 **PRAXISDIENST**
Medical Supplies since 1953









 **Order here!**




vivatmo
me

Instructions for use

English2



Oxycap



Content

1	Your Vivatmo me	4
1.1	Welcome to <i>Vivatmo me</i> , the monitoring system for your respiratory disease	4
1.2	Activate your <i>Vivatmo me</i>	5
1.3	Intended use/ Indications for use	6
1.4	User interface	8
2	Measurement of FeNO in exhaled breath	10
2.1	Attaching the disposable mouthpiece	10
2.2	Performing the measurement	11
2.3	Measurement failed	14
2.4	Interpretation of results	15
2.5	Displaying stored measured values	15
3	Settings/Transmission of Measurement Data	16
3.1	Deactivate/ Activate <i>Bluetooth</i>	16
3.2	Activating/ deactivating audio sound	17
3.3	Changing date and time	18
4	Advanced Settings	19
4.1	Setting the upper and lower thresholds for the traffic light scale	19
4.2	Clearing the memory / reset settings	20

5	Maintenance and Troubleshooting.....	21
5.1	Battery indicator	21
5.2	Cleaning and Disinfecting	21
5.3	Storage	22
5.4	Reconditioning	23
5.5	Quality Control	24
5.6	Troubleshooting	25
6	Technical information.....	27
6.1	Device data	27
6.2	EMC information	30
6.3	Symbols	30
6.4	Disposal of device, batteries and mouthpieces	32
6.5	Warranty	33
7	Manufacturer	33

1 Your Vivatmo me

1.1 Welcome to Vivatmo me, the monitoring system for your respiratory disease

Please consider before using *Vivatmo me*:

- Take the time to read these instructions for use carefully before using the monitoring system.
You need to understand the functions, warnings, displays and operations for a safe and reliable use.
- Set the time and date on your *Vivatmo me* before you begin a measurement.

WARNING

Do not use any defective devices or accessories to avoid malfunction or hazardous situations.


If you need service, information, have a problem or want information about data privacy, please visit the Bosch Vivatmo website at www.vivatmo.com or contact your distributor or the *Vivatmo me* customer service center at service@vivatmo.com. Please have your serial number available when contacting the service center, see 12-digit code on your *Vivatmo me* next to the **SN** symbol at the bottom of the device.

Please keep these instructions for future reference.


Bosch privacy statement:

This product stores measurement data, but no personal information about the user. Configuration, changes and occurring risks by connection to other devices via *Bluetooth®* are in the responsibility of the user of *Vivatmo me*.

1.2 Activate your Vivatmo me



1. Open the battery compartment on the back of the *Vivatmo me* as operator.
2. Insert the 4 AAA primary batteries as shown on the inside of the battery compartment.



3. Remove the battery cover.
4. If the batteries are inserted correctly, the device will automatically turn on.
5. After turning on all display icons will be shown and the system displays the software-version.
6. Symbols to attach disposable mouthpiece are shown.



WARNING

For safe use:

- Turn device off before battery replacement if batteries are inserted.
- Use batteries from the same type as delivered with the device (lithium/iron disulfide or alkaline cells).
- The date and time setting will remain stored in the memory for 5 minutes when the batteries are removed. However, check whether the time and date setting is still correct after the battery replacement.
- Remove the batteries if you are not using the device for a long period of time. The stored measured values are retained.

1.3 Intended use/Indications for use

Intended use: The *Vivatmo me* is intended for quantitative measurement of fractional nitric oxide in human breath.

Indication: Measurement of changes in the fractional nitric oxide concentration in expired breath aids in evaluating a patient's response to anti-inflammatory therapy, as an adjunct to established clinical and laboratory assessments of inflammatory processes such as asthma.

The Bosch *Vivatmo me* system is a non-invasive self-testing device intended to be used at home (in-vitro diagnostic use) as an aid to monitor airway diseases whilst under the care of their physician or healthcare expert. *Vivatmo me* measurement procedure requires patient's cooperation by breathing into the device via a disposable mouthpiece. Patients should be 7 years or older and capable to complete the breathing maneuver.

The measurement procedure of the *Vivatmo me* system generates the fraction of the exhaled breath (FeNO) based on the recommendations for the measurement of exhaled respiratory nitric oxide of the European Respiratory Society (ERS) and the American Thoracic Society (ATS). FeNO is recommended by the ATS in the diagnosis of eosinophilic airway inflammation and in determining the likelihood of responsiveness to anti-inflammatory pharmacological therapy in individuals with chronic respiratory symptoms possibly due to airway inflammation [ATS, 2011]. *Vivatmo me* should only be used as directed in the *Vivatmo me* Instructions for use and as recommended by healthcare professionals. Regardless of displayed test results, if signs or symptoms of chest tightness, shortness of breath, coughing or wheezing occur, the user should contact their healthcare professional immediately.

Contra-indications: Not known.

WARNING

For individual assignment mark device, if several devices are used in the same household. When the device shall be used by another patient, continuously stored data must be deleted to avoid misinterpretation.

The following conditions can influence correct measurement results or cause malfunction or disturbance to the *Vivatmo me* and the disposable mouthpiece:

- Smoking or tobacco consumption for at least 30 minutes before the measurement.
- Eating or drinking shortly before the measurement.
- Strenuous exercise.
- Rooms with high air pollution or open flames, for example gas range or tobacco smoke.
- Ingress of moisture or liquid.

Device contains small parts which can lead to choking.

Children shall use the device supervised by an adult only to ensure correct use of the device.

Store the device outside the access of children when not in use.

When liquid ingresses into the device, remove batteries and prevent further usage of the device to avoid malfunction.

Literature

ATS / ERS, 2005: ATS / ERS recommendations for Measurement of Exhaled Respiratory Nitric Oxide.

ATS, 2011: An Official ATS Clinical Practice Guideline: Interpretation of Exhaled Nitric Oxide Levels (FeNO) for Clinical Application.

1.4 User interface


- ① Traffic light scale – a bullet point displays the result on the scale
Ask your physician or healthcare professional for interpretation advise
- ② Reconditioning
- ③ Trend since last measurement:
 - ↑ The measured value has increased by at least 10 ppb (parts per billion)
 - ↓ The measured value has decreased by 10 ppb (parts per billion) or more
- ④ Invalid disposable mouthpiece: Use a new disposable mouthpiece
- ⑤ Attach disposable mouthpiece
- ⑥ Measuring ID
- ⑦ Battery status
- ⑧ Change date and time or device busy
- ⑨ Activate / deactivate audio sound
- ⑩ Activate / deactivate *Bluetooth* connection
- ⑪ Time and date of the performed measurement
- ⑫ Device ready for measurement
- ⑬ Measured value

- Ⓐ Mouthpiece: Place your lips tightly around the mouthpiece and exhale into the *Vivatmo me*
- Ⓑ Display: Shows current and stored measuring data and device messages
- Ⓒ ON/OFF switch: Activate and deactivate the device
- Ⓓ ENTER button: Press this button to confirm the selection in the display
- Ⓔ ARROW buttons: Press these buttons to change the settings and to access stored measuring data
- Ⓕ Breath intensity display LED: During the measurement, you get a feedback on your breath intensity
 - Blue:** Your breath intensity is right
 - Yellow:** Caution, your breath intensity is too weak
 - Red:** Warning, your breath intensity is too strong
- Ⓖ Battery compartment: Open the battery compartment to exchange the batteries if required

2 Measurement of FeNO in exhaled breath

2.1 Attaching the disposable mouthpiece

The *Vivatmo me* disposable mouthpiece chemically prepares the breath test sample. Human breath is cleaned from contaminants and stabilized to guarantee reliable test results. Use the disposable mouthpiece directly after opening the pouch. The disposable mouthpiece must be attached firmly to perform a measurement.



1. Remove the protective cap from the device.
2. Open the pouch of the disposable mouthpiece.
3. Check the device and disposable mouthpiece and do not use when any damage is noted.



4. Attach the disposable mouthpiece firmly and make sure that it fits securely on the device.
The disposable icon on the display stops moving.



5. "Invalid mouthpiece" is shown on the display if more than 5 attempts have been performed with the same mouthpiece or if the filter capability is used up by the reconditioning.
Remove the mouthpiece and attach a new one.

WARNING

Disposable mouthpiece shall be handled with care for correct measurement results:

- Only use disposable mouthpieces approved for this device.
- Avoid leakage when fitting the disposable to the device.
- Do not clean or dry the disposable mouthpiece.
- Use each disposable mouthpiece for one person and one measurement only, 5 attempts are maximal possible with 1 mouthpiece. Use the mouthpiece within 15 minutes after opening its pouch.
- Take care to use regular disposable mouthpieces with transparent cover. Level-0 disposable mouthpieces with white cover lead to 0 ppb measurements, see chapter "5.5 Quality Control".
- Use only disposable mouthpieces from unopened and undamaged pouches that has not expired. See the expiration date on the package.

2.2 Performing the measurement

For measuring the FeNO value, you have to exhale through the device with a constant breathing rate:

1. Turn ON the device, see "1.2 Activate your *Vivatmo me*".
2. Attach a new disposable mouthpiece, see chapter "2.1 Attaching the disposable mouthpiece".
3. Press the ENTER button to start the regeneration process and wait until the regeneration process is completed. This may take up to 100 seconds. During the regeneration process segments of the display will be animated to indicate the regeneration process is working and the remaining regeneration time is shown.




If the *Vivatmo me* device has been stored in warm and humid environment or not been in use for a longer time, the wrench symbol may appear and a reconditioning will start automatically to provide consistent measurement performance, see chapter “5.4 Reconditioning”.



4. When the display shows the exhaling icon, perform the measurement within 60 seconds (see next steps). The LED on the back of the device turns blue.



5. Turn the device display downwards and make sure you can see the LED at the back of the device. Sit in a relaxed position, inhale deeply through the mouth and lift the device to your mouth. Avoid inhalation through the nose, to exclude influence of nasal NO concentrations.
6. Exhale through the disposable mouthpiece in a controlled manner like blowing into a flute or cooling a hot drink for 10 seconds. Feedback given by the LED on the back of the device:



LED color	Breathing advise for your exhalation
blue	correct
yellow	too weak, caution
red	too strong, warning



If the measurement is successful, the LED flashes blue and the result appears on the display.

The result "5 ppb -LO-" represents a result below the lower detection limit of the device.

The result "300 ppb -HI-" represents a result above the upper detection limit of the device.

7. All results are automatically stored on the device. If the measurement fails, the icon "exhalation failed" appears on the display.

8. Remove and discard the disposable mouthpiece, see chapter "6.4 Disposal of device, batteries and mouthpieces".

9. Place the protective cap on the device.

To save battery power hold the ON / OFF button until the display turns off. The device might show being busy with the clock icon blinking. If you press the ON / OFF button in this mode, the device will turn off afterwards. The device will turn off automatically after 5 minutes of inactivity.

2.3 Measurement failed

Measurement failed when

- the exhalation through the device was excessively weak or strong.
- the measurement was not performed within 60 seconds after completed regeneration process.



Display shows exhalation failed, the LED flashes red and an audio signal occurs.

The regeneration process must be repeated for a new measurement trial:

1. If the clock icon is blinking, the device is busy. Wait until the blinking stops.
2. Press the ENTER button.
3. When the device is ready to measure, perform steps 3 to 7 of the procedure, see chapter "2.2 Performing the measurement".

2.4 Interpretation of results

The measured value is shown on the display. The bullet point under traffic light scale ranks the result relative to the patient's threshold references.

INFORMATION

Interpretation of the results and setting of thresholds must be advised by your physician or healthcare professional.

Default values set for the traffic light scale for adults in accordance with the ATS Guideline 2011, see "Literature" in chapter 1.3:

Green: < 25 ppb, Yellow: 25 ppb - 50 ppb, Red: > 50 ppb

For children below 12 years ATS guideline defines:

Green: < 20 ppb, Yellow: 20 ppb - 35 ppb, Red: > 35 ppb

If you have further questions, contact your distributor or the *Vivatmo me* service center.

2.5 Displaying stored measured values

The *Vivatmo me* automatically stores up to 1,000 measured FeNO values identified by time, date, and ID number. The measured values are stored in descending order starting with the latest. You can retrieve the measured values anytime.

- Press one of the ARROW buttons to show older measured values.
- Press the ENTER button to exit the view.

3 Settings / Transmission of Measurement Data

3.1 Deactivate / Activate *Bluetooth*

Measurement data can be transmitted encryptedly via *Bluetooth* to a remote device like a smartphone with the *Vivatmo app* (see www.vivatmo.com).



1. Change *Bluetooth* by holding the ENTER button for 2 seconds.
The settings toolbar appears.



2. In the menu screen, use the ARROW buttons to select the *Bluetooth* icon.
The icon starts flashing.

3. Press the ENTER button to activate / deactivate *Bluetooth*.
Press one of the ARROW button to confirm selection.

4. To exit the basic settings toolbar press the left or right ARROW button repeatedly until icons stop flashing, then press the ENTER button.

3.2 Activating/deactivating audio sound



1. Hold the ENTER button for 2 seconds. The settings toolbar appears.



2. In the menu screen, use the ARROW buttons to select the Audio icon. The selected icon flashes.

3. Press the ENTER button to activate / deactivate the sound. If the check mark is shown, the audio is on. Press one of the ARROW buttons to confirm selection.

4. To exit the basic settings toolbar press the left or right ARROW button repeatedly until no icons are flashing, then press the ENTER button.


3.3 Changing date and time



1. Hold the ENTER button for 2 seconds.
The settings toolbar appears.




2. In the menu screen, use the ARROW buttons to select the Clock icon, the icon starts flashing, then press the ENTER button.



3. Use the ARROW buttons to choose the date format dd.mm.yy or mm/dd/yy, then press the ENTER button.



4. Use the ARROW buttons to select the year, then press the ENTER button.
5. Use the ARROW buttons to select the month, then press the ENTER button.
6. Use the ARROW buttons to select the day, then press the ENTER button.



7. Use the ARROW buttons to switch between 0 - 24 h and 0 - 12 h, AM / PM time format, then press the ENTER button.



8. If you have chosen 0 - 12 h time format, you can select between AM and PM, then press the ENTER button.
9. Use the ARROW buttons to set the hours, then press the ENTER button. Then set the minutes and press the ENTER button.
10. To exit the basic settings toolbar, press the left or right ARROW button repeatedly until icons stop flashing, then press the ENTER button.

4 Advanced Settings

WARNING

These settings shall be performed by healthcare professionals only. Incorrect operation may lead to misinterpretation of the measuring results or deletion of the stored data.
If you have further questions, contact the *Vivatmo me* service center.



In order to access the advanced settings, press the ENTER button and the right ARROW button for 3 seconds.

4.1 Setting the upper and lower thresholds for the traffic light scale

The default values set for the traffic light scale are the reference values for adults in accordance with the ATS Guideline 2011, see chapter “2.4 Interpretation of results”.



1. Use the ARROW buttons to select the Thresholds icon. The icon starts flashing.
Press the ENTER button to select the flashing icon
2. Use the ARROW buttons to set the threshold value for the lower threshold. Confirm with ENTER button.



3. Use the ARROW buttons to set the threshold value for the upper threshold. Confirm with ENTER button.
4. To exit the advanced settings toolbar press the ARROW buttons until no icons are flashing, then press the ENTER button to return to the main screen.



4.2 Clearing the memory/reset settings



1. Use the ARROW buttons to select the Clear memory icon. The icon start flashing.
Press the ENTER button to select the flashing icon.
2. Use the ARROW buttons to select the check mark.
3. Press the ENTER button to confirm the clearing.
The Clear memory icon is shown and flashes.
4. To exit the advanced settings toolbar, press the ARROW buttons until no icons are flashing, then press the ENTER button to return to the main screen.

5 Maintenance and Troubleshooting

5.1 Battery indicator

no icon		The batteries charge is sufficient.
	shown	The charge is low and the batteries should be replaced soon.
	blinking	The batteries should be replaced before using the device again.

INFORMATION

When replacing the batteries, always replace the whole set of batteries, see chapter "1.2 Activate your *Vivatmo me*".

5.2 Cleaning and Disinfecting

Clean the device regularly or when used by another patient.

1. Turn off the device. Wash hands thoroughly with soap and water.
2. Wipe the entire *Vivatmo me* with disinfectant wipes with maximum 30 % alcohol e.g. Mikrobac Tissues (Bode Chemie GmbH) or Sagrotan wipes (Reckitt Benckiser) or with a soft cloth, moistened with soap cleaning agents. Carefully wipe around openings.
3. For disinfecting repeat step 2 with disinfectant wipes with maximum 30 % alcohol following the wipe's manufacturer's instructions for disinfecting.

WARNING

The following conditions can cause malfunction to the *Vivatmo me* or the disposable mouthpiece:

- Fractures in the housing.
- Ingress of spray, moisture or liquid.

5.3 Storage

Protect your *Vivatmo me*:

- Always attach the protective cap when *Vivatmo me* is not in use.
- Keep the *Vivatmo me* and its components in a clean, cool and dry place.

WARNING

The following conditions can cause damage and may cause malfunction to the *Vivatmo me* and the disposable mouthpiece:

- Storage in high humidity, high temperature, direct sun radiation or high exposure to dust or volatile substances (e. g. disinfectants or nail polish remover).
- Places that are subject to vibration, shock, or near hot surfaces.
- Rooms with open flames (e. g. gas range) or smoke (e. g. from tobacco consumption).

5.4 Reconditioning

If the *Vivatmo me* device has been stored in warm and humid environment or not been in use for a longer time, the wrench symbol may appear and a reconditioning is required to provide consistent measurement performance. In this case the wrench symbol appears during the regeneration process on the pacifier screen and the device automatically starts the reconditioning.



1. The wrench symbol indicates that the reconditioning process is working.
Segments of the display will be animated and the remaining reconditioning time is shown in minutes.



2. The system beeps and blinks 3 times blue when reconditioning is finished and the device shows the Exhaling icon when it is ready to measure. Perform steps 3 to 7 of the procedure, see chapter "2.2 Performing the measurement".
3. Use the ARROW buttons to view the stored data during the reconditioning and return with ENTER to the animated screen.

INFORMATION

The reconditioning takes between 7 up to 90 minutes and cannot be interrupted. The wrench symbol is shown during all activities:

- Change the batteries if E-6 and the blinking battery symbol appears during reconditioning, see chapter "1.2 Activate your *Vivatmo me*".
- Change the disposable mouthpiece if the "Invalid mouthpiece" screen appears during reconditioning, see chapter "2.1 Attaching the disposable mouthpiece".
- If you push the ON/OFF button, OFF appears as long as the reconditioning is active.

5.5 Quality Control

Quality Control lets you know that your device gives reliable results.

Perform quality control when:


- You dropped the device.
- You want to make sure you are testing correctly.

Confirm a measurement below the detection limit of 5 ppb by using a white Level-0 disposable mouthpiece. Above the detection limit compare with results achieved with another certified FeNO measurement system under consideration of the measurement performances. When used in clinical settings perform quality control regularly. If one of the confirmations failed, prevent further usage of the device to avoid misinterpretation.

5.6 Troubleshooting

When using the *Vivatmo me*, an error message may indicate a problem.

The following table gives explanations of the error messages as well as recommendations for solving the problem.

Display	Description	Solution
	Regeneration or reconditioning process is ongoing.	Wait for the process to be completed. Do not remove the mouthpiece.
E-1	Permissible number of measurements < 20.	See display for number of remaining measurements.
E-2	Maximum number of permissible measurements reached. Device does not perform anymore measurements.	Dispose <i>Vivatmo me</i> , see chapter "6.4 Disposal of device, batteries and mouthpieces".
E-4	Humidity too high.	Check ambient conditions and take a new mouthpiece, see chapter "2.1 Attaching the disposable mouthpiece".
E-5	Ambient temperature or air pressure is outside the specified range.	Use the <i>Vivatmo me</i> in location with permissible ambient conditions. When transporting the device from one place to another, a stabilization period of min. 20 minutes should be allowed before use.

Display	Description	Solution
E-6	Batteries are too weak for measurement.	Replace batteries with new ones, see chapter "1.2 Activate your <i>Vivatmo me</i> ".
E-7	System self-test failed.	Replace batteries with new ones, see chapter "1.2 Activate your <i>Vivatmo me</i> " and start again. If the error persists, contact the service.
E-8	Regeneration failed.	Press the ENTER button to repeat the regeneration. If the error persists, contact the service.
None – The device will not turn on.	Batteries are empty.	Replace batteries with new ones, see chapter "1.2 Activate your <i>Vivatmo me</i> ".
	Batteries are inserted the wrong way.	Ensure the batteries are inserted with correct polarity.

If you have any problems which you cannot solve with these Instructions for use or the device shows unexpected operations, contact your distributor or the *Vivatmo me* customer service center at service@vivatmo.com. For the contact data see chapter "7 Manufacturer". Please indicate the following information:

- Your name, address and phone number
- Serial number of the device (12-digit code on your *Vivatmo me* next to the **SN** symbol at the bottom of the device)
- Error message
- A detailed description of the problem

6 Technical information

6.1 Device data

Product description	FeNO monitoring system
Model	<i>Vivatmo me</i> , catalogue number F09G100149
Measuring range	5 ppb to 300 ppb
Accuracy	± 5 ppb below 50 ppb, ± 10 % for ≥ 50 ppb, ± 15 % for ≥ 160 ppb expressed as the upper/lower confidence limit of 95 %
Precision	± 5 ppb below 50 ppb, ± 10 % for ≥ 50 ppb, ± 15 % for ≥ 160 ppb expressed as one standard deviation for replicate measurements with the same instrument
Lifetime, device service life, memory capacity	1,000 measurements
Power source	4 AAA batteries 1.5 V, Useful life of primary batteries: Lithium/iron disulfide up to 60 attempts, alkaline up to 25 attempts, (Battery and compartment lid are detachable parts)

Applied part	Type BF as per EN 60601-1 for handheld and disposable mouthpiece when attached
Electrical safety	ME device with internal supply, tested as per EN 60601-1-11 IP22 (protection against solid particles $\geq 12,5$ mm and ingress of dripping water when tilted up to 15°) for basic safety but not for function
Maximum surface temperature	58 °C, touch time < 60 seconds
Electromagnetic emissions	CISPR 11 Group 1 (battery operated)
Electromagnetic immunity	EN 61000-4-2, EN 61000-4-3 (battery operated), EN 61000-4-8
Sensor	Chemical field-effect transistor
Data transfer	<i>Bluetooth</i> Smart (low energy), 2.4 GHz frequency band

Operating temperature / humidity / air pressure	+15 °C to +27 °C / 15 % to 60 % relative humidity (non-condensing) / 780 hPa to 1100 hPa, corresponds to 0 to 2,000 m above sea level
Storage and transport temperature / humidity / air pressure	+5 °C to +27 °C / 10 % to 60 % relative humidity between uses / 780 hPa to 1100 hPa, corresponds to 0 to 2,000 m above sea level
Weight / Dimensions	170 g / 4.0 cm x 5.4 cm x 22.4 cm
Package contents	<i>Vivatmo me</i> device, 5 disposable mouthpieces (Oxycaps), 4 batteries, protective cap (detachable part), Instructions for use
Disposable mouthpiece single use (accessory)	Measurement within 15 minutes after opening the pouch. Use is limited to 5 attempts and expiration date. Lifetime disposable 2 years from date of manufacturing. Further purchase: <i>Vivatmo Oxycap</i> mouthpieces (catalogue number F09G100152), Level 0 mouthpieces (catalogue number F09G100174).

Limitations of the System: Exchange your *Vivatmo me* at the latest 3 years after manufacturing date.
Information concerning candidates of the REACH regulation can be found at www.vivatmo.com.

6.2 EMC information



Important information regarding electromagnetic compatibility (EMC)



Hereby, Bosch Healthcare Solutions GmbH, declares that the radio equipment type *Vivatmo me* is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address: <https://www.vivatmo.com/en/service-support/additional-information.html>






This device complies with EN60601-1-2:2015 for EMC with the objective to avoid insecure product situations. This standard regulates the levels of immunity against electromagnetic interferences and the maximum electromagnetic emission values for medical equipment. This medical device manufactured by the company complies with the standard EN60601-1-2:2015 both in terms of immunity and of emissions and does therefore not need any service and maintenance regarding EMC and ESD over lifetime. For manufacturers declaration regarding EMC see www.vivatmo.com.






Please note that portable and mobile HF communication systems may interfere with this device even if compliant with CISPR emission requirements. Do not staple the device or use any mobile phones or other devices generating strong electrical or electromagnetic fields. This could result in malfunction of the medical device and may create a potentially insecure situation. Portable RF communication devices are not to be used closer than 30 cm next to the device.

6.3 Symbols

	Temperature range
	Application range humidity

	Application range air pressure
	Serial number

	Manufacturer address
	The IVD product meets the requirements of applicable European directives.
	Consult instructions for use
WARNING	Important information to avoid a hazardous situations
INFORMATION	Important advise
 YYYY-MM-DD	Use-by date
	For single-use only
LOT	Batch number

REF	Reference number / type part number
	Applied part type BF as per EN60601-1
IVD	In vitro diagnostic medical device 98 / 79 / EEC IVD Directive
	Keep dry
IPN₁N₂	IP Protection class
	<i>Bluetooth</i> compatible
	Manufacturing date
	Regulatory compliance mark Australia

The *Bluetooth*® word mark and logos are registered trademarks owned by *Bluetooth* SIG, Inc. and any use of such marks by Bosch Healthcare Solutions GmbH is under license.

6.4 Disposal of device, batteries and mouthpieces

INFORMATION

Delete the data on the device prior to disposal (see chapter 4.2).



For the purpose of disposal, please separate this device from other waste to prevent possible harm to the environment or human health from uncontrolled waste disposal. Turn in the device or recycle in order to promote the sustainable reuse of material resources (in compliance with EU directive 2012 / 19 / EU). The used or expired device and batteries should be recycled in compliance with the local recycling program for electronic equipment. Dispose used batteries in compliance with the national / local regulations for the disposal of batteries separated from the measuring unit (in compliance with the European Guideline 2006 / 66 / EC).

The mouthpiece is a single-use product. It contains small quantities of potassium permanganate (KnMnO_4). Use a new mouthpiece for each measurement and dispose the used one in residual waste.

6.5 Warranty

The statutory provisions on warranty rights in consumer goods sales in the country of purchase shall apply. Batteries are not covered by warranty.

WARNING

Do not disassemble *Vivatmo me* or its components, and do not try to repair it by yourself. Any changes or modifications to *Vivatmo me* that are not approved by the manufacturer may cause malfunction and result in the warranty becoming invalid.

7 Manufacturer

Bosch Healthcare Solutions GmbH
Stuttgarter Str. 130
71332 Waiblingen, Germany

E-Mail: service@vivatmo.com
Internet: www.vivatmo.com





ผู้ใดมีไว้ในครอบครอง หรือ ใช้ซึ่งเครื่องวิทยุ
คมนาคมหรือตั้งสถานีวิทยุคมนาคมนี้ ต้องได้รับ
ใบอนุญาตจากเจ้าพนักงานผู้ออกใบอนุญาต
หากฝ่าฝืน มาตรา 6 หรือมาตรา 11 มีความผิด
ตามมาตรา 23 แห่งพระราชบัญญัติวิทยุ
คมนาคม พ.ศ. 2498 ต้องระวางโทษปรับไม่เกิน
หนึ่งแสนบาท หรือจำคุกไม่เกินห้าปี หรือทั้งปรับ
ทั้งจำ



nanb. | โทรคมนาคม
กำกับดูแลเพื่อประชาชน
Call Center 1200 (InswS)