



Technical data sheet

Item No.:	130666
Product name:	Underpad
Product description:	Underpad for incontinence usage, mostly in hospital, senior homes, homecare and so on.
Sizes of pad:	40 x 60 cm
Inner absorbent pad size:	40 x 50 cm
Raw Materials:	Nonwoven, Tissue, PE, Fluff pulp, SAP
Total weight per piece:	40 g (14.35 g Amount of top nonwoven+tissue+ PE) (23.65 g Amount of fluff pulp) (2 g Amount of SAP) (22 gms Back layer PE) (13 gms Tissue) (12 gms Top layer nonwoven)
Absorption:	0.5 liters absorption capacity
Color:	Green
Packaging:	10 pcs / bag 10 bags / carton
Packing Materials:	Polyethylene 40 µm thickness
Carton size:	34 x 29.5 x 30 cm
Useful Life:	5 years
Handling and Storage:	Store dry, Avoid direct sunlight
Storage temperature:	temperature between -5-40 °C

Cautions, Warnings or contraindications

No. 70 Yong'an Road, Xuguan Industrial Park, High-Tech District, Suzhou, Jiangsu, China 215151
Telephone: 86-512-66653331, Fax: 86-512-66651066
www.sn-anpa.com



SUZHOU SUNING UNDERPAD CO., LTD
蘇州市蘇寧床墊有限公司

Not to use if the product is damaged or dirty, or do have other functional failures.

Rubbish and Final Disposition

The product must be disposed, disposable waste regulation, after using.



Date: 2021-11-12