

RESCUE SPACE BLANKET



DESCRIPTION

This rescue blanket has a silver and a gold side to protect against cold or heat. To protect from the heat, the silver side should be facing outward, while the blanket is loosely covering the patient. To protect the patient from the heat, the silver side should be facing outward and the loosely covering the patient. The small size and light weight make this rescue blanket ideal for any emergency bag or first aid kit.

PRODUCT DETAILS

- Made of temperature tolerant foil
- Usable as heat or cold protection
- Reflects body warmth
- Dimensions: 210 cm x 160 cm
- Dimensions (folded): 12.5 cm x 8 cm x 1.8 cm
- Colour: silver/gold
- Delivery contents: 1 piece

Overview

Article Number:
T138113



Category:
Emergency Equipment

EAN:
4260626262060



Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

