



**Choose the right stethoscope to help you perform at your best.**

3M™ Littmann® Stethoscopes feature adjustable headsets with soft-sealing eartips. Select Littmann stethoscopes are designed with tunable technology to give you the ability to hear different frequency sounds by simply adjusting the pressure applied to the chestpiece. Get the reliable acoustic performance, high quality and comfort you need with a Littmann stethoscope.



	Lightweight II S.E. Stethoscope	Classic II Pediatric and Infant Stethoscopes	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	CORE Digital Stethoscope
	Ideal for taking blood pressure readings and making limited patient assessments.	Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.		Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments.		High-performance auscultation. Hear more with amplification, noise cancellation and mobile connectivity.
<b>Approximate Length</b>	28"/71cm	28"/71cm	27"/69cm	22"/56cm, 27"/69cm	27"/68cm	27"/69cm
<b>Approximate Weight</b>	118g	105g, 95g	150g	167g, 177g	185g	232g
<b>Warranty</b>	2 years	3 years	5 years	7 years	7 years	2 years
<b>Diaphragm</b>	Tunable	Non-Tunable (Floating)	Tunable (2)	Tunable (2)	Tunable	Tunable (2)
<b>Approximate Chestpiece Diameter</b>	Double-sided 2.0 x 1.8"/5.2 x 4.5cm 1.4 x 1.3"/3.5 x 3.2cm	Double-sided Pediatric 1.4"/3.4cm 1.1"/2.7cm Infant 1.1"/2.8cm 0.7"/1.7cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm	Single-sided 1.9"/4.8cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm
<b>Features</b>	Teardrop-shaped chestpiece Fits under blood pressure cuff	Specially designed for the smallest patients	Convertible open bell	Convertible open bell Dual-lumen tubing	Dual-lumen tubing Best acoustic response in Littmann mechanical line	Up to 40x amplification* Switch from analog listening to digital  Connects to Eko software on a smart device** to capture, store and share data



**The power to learn, in the palm of your hand.**

3M™ Littmann® Learning App – Download free version in the App Store® or Google Play™ Store

**Need more help? Visit [Littmann.com.au/FindMyStethoscope](http://Littmann.com.au/FindMyStethoscope)**  
**Need more help? Visit [Littmann.co.nz/FindMyStethoscope](http://Littmann.co.nz/FindMyStethoscope)**

1300 363 878 | [www.littmann.com.au](http://www.littmann.com.au)  
0 800 80 8182 | [www.littmann.co.nz](http://www.littmann.co.nz)

\*At peak frequency, vs. analog mode  
\*\*Not included



# 3M™ Littmann® Stethoscope Ordering Information and Tube Colour Guide

**Components**  
(Chestpiece/Eartube/Stem)

Lightweight II SE Stethoscope	Classic II Infant Stethoscope	Classic II Pediatric Stethoscope	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	CORE Digital Stethoscope
-------------------------------	-------------------------------	----------------------------------	-------------------------------------	---------------------------------------	--------------------------------	--------------------------

**Tube Colour**

All Standard Finish	2450 Black	2114 Black	2113 Black	5620 Black	5627 Burgundy	6156 Plum	2160 Black
	2451 Burgundy	2114R Red	2113R Red	5621 Gray	5633 Pearl Pink	6152 Black	2163 Burgundy
	2452 Caribbean	2124 Caribbean	2119 Caribbean	5622 Navy	5831 Plum	6154 Navy	2164 Navy
	2454 Ceil Blue		2122 Raspberry	5623 Caribbean	5832 Lavender	6155 Hunter Green	2167 Plum
	2456 Pearl Pink			5630 Ceil Blue	5835 Turquoise		
				5839 Lemon-Lime			
All Black			5803 Black	5868 Burgundy	6163 Black	6168 Navy	2161 Black
All Brass Finish							8480 Black
All Copper Finish			5809 Chocolate				
All Smoke			5811 Black				
Rainbow/Brass/Brass		2157 Raspberry		5806 Raspberry	5807 Caribbean		
Brass/Smoke/Smoke							
Rainbow/Black/Black			5870 Black				
Mirror/Black/Black							8890 Black
Mirror/Smoke/Smoke				5863 Navy			
			5959 Ceil Blue				
Mirror/Standard/Mirror					6177 Black	6170 Burgundy	
Champagne/Smoke/Smoke			5861 Black	5864 Burgundy	6179 Black		
Champagne/Black/Champagne					6176 Burgundy	6190 Caribbean	
High-Polish Copper/Black/Black							8870 Black
High-Polish Rainbow/Black/Black					6242 Navy		8570 Black
High-Polish Rainbow/Smoke/Smoke					6240 Black	6241 Raspberry	
<b>COLOURED STEM SERIES</b>	<b>COLOURED STEM SERIES: (square inset represents stem colour)</b>						
Black/Black/Colour	<b>Components:</b> Chestpiece/Eartube/Stem					6200 Black	6203 Black
High-Polish Smoke/Black/Colour	<b>Standard Finish</b> 	<b>Rainbow</b> 				6202 Navy	6204 Black
High-Polish Smoke/Smoke/Colour	<b>Black</b> 	<b>Mirror</b> 				6232 Black	6234 Caribbean
High-Polish Rainbow/Black/Colour	<b>Brass Finish</b> 	<b>Champagne</b> 				6238 Gray	
Rainbow/Black/Colour	<b>Copper Finish</b> 	<b>High-Polish Smoke</b> 					
Smoke/Smoke/Colour	<b>Smoke</b> 	<b>High-Polish Rainbow</b> 					
Mirror/Smoke/Colour		<b>High-Polish Copper</b> 		5872 Turquoise	5875 Lime		
				5873 Gray			
				5960 Plum	5962 Pearl Pink		

Colours shown are approximate.

1300 363 878 | [www.littmann.com.au](http://www.littmann.com.au)  
0800 80 81 82 | [www.littmann.co.nz](http://www.littmann.co.nz)

© 3M 2021. All rights reserved. 3M, Littmann, the L Littmann logo, Classic III, Cardiology IV, and Master Cardiology are trademarks of 3M. Eko trademarks, including CORE, are marks and/or registered marks of Eko Devices Inc. All other trademarks are property of their respective owners. 70-2011-8157-8, Rev. 28/01/22