


COVID-19 Antigen Rapid Test Cassette



Nasal Swab



(For Self-Testing)

The COVID-19 Antigen Rapid Test Cassette(Nasal Swab) is a rapid chromatographic immunoassay for the qualitative detection of COVID-19 antigen in Nasal Swab in symptomatic individuals. The identification is based on the monoclonal antibodies specific for the Nucleocapsid (N) protein of SARS-CoV-2. It is intended to aid in the rapid differential diagnosis of current COVID-19 infections.

Painless

— Avoid the discomfort of nasopharyngeal swab sampling

Convenient

— Easy to operate, Interpret result at 10 minutes

Accurate

— Excellent performance compared to molecular methods

 **Order here!**

 **PRAXISDIENST**
Medical Supplies since 1953

Hangzhou Biotest Biotech Co., Ltd.


17#, Futai Road, Zhongtai Town, Yuhang District, Hangzhou -311121 P.R. China
TEL: +86-571-89058016 • FAX: +86-571-88633388 • EMAIL: info@biotests.com.cn • Website: en.biotests.com.cn

COVID-19 Antigen Rapid Test Cassette

Materials

- Test cassette
- Extraction Buffer Tube
- Sterile Nasal Swab
- Package Insert
- Quick Reference Guide
- Disposal bag
- Workstation (25T Only)

Convenient Procedure



Ordering Information



Cat.No.	Product	Specimen	Format	Pack
ICOVN-C81H	COVID-19 Antigen Rapid Test Cassette (For Self-testing) ^{CE1434}	Nasal Swab	Cassette	1 T (soft packaging)
ICOVN-C81H	COVID-19 Antigen Rapid Test Cassette (For Self-testing) ^{CE1434}	Nasal Swab	Cassette	1 T
ICOVN-C81H	COVID-19 Antigen Rapid Test Cassette (For Self-testing) ^{CE1434}	Nasal Swab	Cassette	2 T
ICOVN-C81H	COVID-19 Antigen Rapid Test Cassette (For Self-testing) ^{CE1434}	Nasal Swab	Cassette	5 T
ICOVN-C81H	COVID-19 Antigen Rapid Test Cassette (For Self-testing) ^{CE1434}	Nasal Swab	Cassette	25 T

Hangzhou Biotest Biotech Co., Ltd.

17#, Futai Road, Zhongtai Town, Yuhang District, Hangzhou -311121 P.R. China
 TEL: +86-571-89058016 • FAX: +86-571-88633388 • EMAIL: info@biotests.com.cn • Website: en.biotests.com.cn