

# Economy anaesthetic face masks

Seven sizes from neonatal to extra large adult



Now  
colour-  
coded!


# Economy anaesthetic face masks

## Introduction

A range of colour-coded economy anaesthetic face masks for single patient use to prevent the possibility of cross infection.

## Features

- Soft cushion for patient comfort and air tight fit.
- Clear flexible housing and cushion for good patient visibility.
- Seven sizes available, from neonatal to extra large adult
- Each mask supplied with colour-coded hook ring.



## Economy anaesthetic face masks

|             |                             |        |             |    |     |             |  |     |
|-------------|-----------------------------|--------|-------------|----|-----|-------------|--|-----|
| <b>1517</b> |                             | 22F    | <b>1516</b> |    | 22F | <b>1515</b> |  | 22F |
| <b>1514</b> |                             | 22F    | <b>1513</b> |    | 22F | <b>1512</b> |  | 15M |
| <b>1511</b> |                             | 15M    |             |    |     |             |  |     |
| <b>1517</b> | Extra large adult face mask | size 6 |             | 25 |     |             |  |     |
| <b>1516</b> | Large adult face mask       | size 5 |             | 30 |     |             |  |     |
| <b>1515</b> | Medium adult face mask      | size 4 |             | 40 |     |             |  |     |
| <b>1514</b> | Small adult face mask       | size 3 |             | 50 |     |             |  |     |
| <b>1513</b> | Paediatric face mask        | size 2 |             | 35 |     |             |  |     |
| <b>1512</b> | Infant face mask            | size 1 |             | 32 |     |             |  |     |
| <b>1511</b> | Neonatal face mask          | size 0 |             | 40 |     |             |  |     |

IS5.6 Issue 5 AAC 01.07

**INTERSURGICAL** COMPLETE RESPIRATORY SYSTEMS United Kingdom (Head Office): Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ  
 T: +44 (0)118 9656 300 F: +44 (0)118 965 6356 info@intersurgical.com | www.intersurgical.com

France  
info@intersurgical.fr

Deutschland  
info@intersurgical.de

España  
info@intersurgical-es.com

Россия  
info@intersurgical.ru

South Africa  
info@intersurgical.co.za


Portugal  
info@intersurgical.pt

Nederland  
info@intersurgical.nl

Lietuva  
info@intersurgical.lt

Česká Republika  
info@intersurgical.cz

Philippines  
info@intersurgical.ph



Certified to ISO 9001 and ISO 14001